

## ***MRI Liver / Biliary / Pancreas - Prohance***

Application Examples: abdomen pain, hemangioma, metastatic disease from primary other than liver

FoV Coverage: Top of diaphragm to inferior liver edge

| Sequence                         | FoV     | Thickness | Gap    |                         |
|----------------------------------|---------|-----------|--------|-------------------------|
| Cor T2 Haste/FSFSE               | 42.0 cm | 6.0 mm    | 0.5mm  |                         |
| Axial T2 FSFSE FS                | 42.0 cm | 6.0 mm    | 0.5mm  |                         |
| Axial DWI                        | 40.0 cm | 8.0 mm    | No Gap | B Values 0,600,800,1000 |
| Axial In/Out Phase               | 40.0 cm | 4.8 mm    | No Gap |                         |
| Axial T1 FS                      | 40.0 cm | 6.0 mm    | 0.5mm  |                         |
| Axial Vibe/Lava Flex             | 40.0 cm | 4.4 mm    | No Gap |                         |
| INJECT CONTRAST (1mmol/Kg)       |         |           |        |                         |
| Axial T1 Vibe/Lava Flex Arterial | 40.0 cm | 4.4 mm    | No Gap | 20 second delay         |
| Axial T1 Vibe/Lava Flex Venous   | 40.0 cm | 4.4 mm    | No Gap | 50 second delay         |
| Coronal Vibe/Lava Flex Venous    | 40.0 cm | 4.4 mm    | No Gap |                         |
| Axial T1 Vibe/Lava Flex          | 40.0 cm | 4.4 mm    | No Gap | 5-minute delay          |