

MRI Liver / Biliary / Pancreas - Eovist

Application Examples: hepatocellular carcinoma in cirrhotic patient, liver lesion, FNH vs adenoma mets specifically to liver
FoV Coverage: Top of diaphragm to inferior liver edge

Sequence	FoV	Thickness	Gap	
Cor T2 Haste/FSFSE	42.0 cm	6.0 mm	0.5 mm	
Axial T2 FSFSE FS	40.0 cm	6.0 mm	0.5 mm	
Axial DWI	38.0 cm	8.5 mm	No Gap	B Values 0,600,800,1000
Axial In/Out Phase	38.0 cm	4.8 mm	No Gap	
Axial T1 FS	38.0 cm	4.8 mm	No Gap	
Axial Vibe/Lava Flex	40.0 cm	4.4 mm	No Gap	
INJECT CONTRAST (1mmol/Kg)-Rate 2cc/sec				
Axial T1 Vibe/Lava Flex Arterial	40.0 cm	4.4 mm	No Gap	20 second delay
Axial T1 Vibe/Lava Flex Venous	40.0 cm	4.4 mm	No Gap	50 second delay
Coronal Vibe/Lava Flex Venous	40.0 cm	4.4 mm	No Gap	120 second delay
Axial T1 Vibe/Lava Flex	40.0 cm	4.4 mm	No Gap	5-minute delay
Sagittal FSFSE/HASTE	42.0 cm	6.0 mm	0.5 mm	
Axial T1 Vibe/Lava Flex	40.0 cm	4.4 mm	No Gap	20-minute delay