Gundersen Health System BODY MR Protocol GE 1.5T 450W

MRI MRCP

Application Examples: elevated ALT, common bile duct dilation, common bile duct stone FoV Coverage: Liver dome to approximately mid kidney to include all of pancreas

Sequence	FoV	Thickness	Gap	
Coronal Haste/FSFSE	35.0 cm	6 mm	0.5 mm	Through abdomen
Axial Haste/FSFSE	35.0 cm	6 mm	0.5 mm	Through abdomen
Sagittal Haste/FSFSE	35.0 cm	6 mm	0.5 mm	Through abdomen
Coronal Haste/FSFSE	35.0 cm	3 mm	1.0 mm	Through abdomen
Axial Haste/FSFSE	35.0 cm	3 mm	1.0 mm	Through abdomen
Sagittal Haste/FSFSE	35.0 cm	3 mm	1.0 mm	Through abdomen
Coronal Thick slab RARE	30.0 cm	40.0 cm	0.0 mm	Through common bile duct
Sagittal Thick slab RARE	30.0 cm	40.0 cm	0.0 mm	Through common bile duct
				Through common bile duct
Coronal 3D RARE	35.0 cm	1.4 mm	0.0 mm	MIP and send source images to PACS