

MRI Abdomen/Pelvis

Application Examples: Mets

FoV Coverage: Top of diaphragm to pubic symphysis

Sequence	FoV	Thickness	Gap	
Coronal Haste/SSFSE	44.0 cm	6.0 mm	0.5 mm	AP
Axial In/Out Phase	40.0 cm	4.4 mm	No Gap	Abdomen only
Axial T2 FS Haste/SSFSE	44.0 cm	6.0 mm	0.5 mm	AP
Axial Vibe FS/LAVA Flex Pre	40.0 cm	4.4 mm	No Gap	AP
Axial DWI Abd/Pelvis	40.0	6.0mm	0.5mm	B Value 0,600,800,1000
INJECT CONTRAST (1mmol/Kg)				
Axial Vibe FS/LAVA Flex post	40.0 cm	4.4 mm	No Gap	20 sec delay Abdomen
Axial Vibe FS/LAVA Flex post	40.0 cm	4.4 mm	No Gap	50 sec delay Abdomen
Axial Vibe FS/LAVA Flex post	40.0 cm	4.4 mm	No Gap	120 second delay AP