

MRI Adrenal w/o contrast

Application Examples: adrenal lesion

Sequence	FoV	Thickness	Gap	Gap
Coronal Haste/SSFSE	44.0cm	6.0mm	0.5mm	Liver dome to mid kidney
Axial T2 FS	38.0cm	6.0mm	0.5mm	Respiratory triggered
Axial In/Out phase	38.0cm	4.8mm	0.5mm	Liver dome to mid kidney
Axial T1 FS	38.0cm	6.0mm	0.5mm	Liver dome to mid kidney