

MRI Uterus or Cervical Cancer 3T only

Application Examples: Cervical cancer, uterine cancer, vaginal cancer

Patient Prep: KY gel in vagina

Sequence	FoV	Thickness	Gap	
Axial T2	26.0 cm	5.0 mm	0.8 mm	Bifurcation to Pubic Symphysis
Axial T1	26.0 cm	5.0 mm	0.8 mm	Bifurcation to Pubic Symphysis
Sagittal T2	18.0 cm	3.0 mm	0.5 mm	From top of Uterus to Pubic Symphysis
Axial T2	18.0 cm	3.0 mm	1.0 mm	Short axis of Uterus/CX only
Coronal T2	18.0 cm	3.0 mm	1.0 mm	Long axis of Uterus/CX only
Axial Diffusion	28.0 cm	4.0 mm	No Gap	Short axis of Uterus/CX only B Values 0, 600, 800, 1000
Sagittal Diffusion	28.0 cm	4.0 mm	No Gap	B values 0, 400,800
Sagittal Vibe FS/Lava Flex Pre	40.0 cm	4.4 mm	No Gap	From top of Uterus to Pubic Symphysis
INJECT CONTRAST (1mmol/Kg)				
Sagittal Vibe FS/Lava Flex	40.0 cm	4.4 mm	No Gap	40 seconds after injection
Sagittal Vibe FS/Lava Flex	40.0 cm	4.4 mm	No Gap	90 Seconds after injection
Axial T1 FS with gad	26.0 cm	5.0 mm	0.8 mm	Bifurcation to Pubic Symphysis