

MRI Uterine Fibroid 3T Only

Application Examples: uterine fibroid

Patient Prep: KY gel in vagina

Sequence	FoV	Thickness	Gap	
Axial T2	36.0cm	5.0mm	0.8mm	Bifurcation to Pubic Symphysis
Axial T1	36.0cm	5.0mm	0.8mm	Bifurcation to Pubic Symphysis
Sagittal T2 FS	25.0cm	5.0mm	1.5mm	From top of Uterus to Pubic Symphysis
Axial GRE (Medic/Merge)	36.0cm	5.0mm	0.8mm	From top of Uterus to Pubic Symphysis
Axial T2 HIREs	18.0cm	3.0mm	1.0mm	Short axis of Uterus/CX only
Coronal T2 HIREs	18.0cm	3.0mm	1.0mm	Long axis of Uterus/CX only
Sagittal Diffusion	25.0	5.0mm	1.5mm	B values 0, 1000
Sagittal Vibe FS/Lava Flex Pre	30.0cm	4.0mm	0.0 mm	From top of Uterus to Pubic Symphysis
INJECT CONTRAST (1mmol/Kg)				
MRA Cor Tricks/Twist	30.0	2.0mm	0.0 mm	Arterial
Sagittal Vibe FS/Lava Flex	30.0cm	4.0mm	0.0 mm	90 seconds after injection
Axial T1 FS with gad	36.0cm	5.0mm	0.8mm	Bifurcation to Pubic Symphysis