

***MRI Uterine Anomaly 3T only***

Application Examples: Uterine anomaly, infertility

Patient Prep: KY gel in vagina

FoV Coverage: Top of uterus to pubic symphysis

Sequence	FoV	Thickness	Gap	
Coronal Haste/FSFSE	45.0 cm	6.0 mm	0.5 mm	Abdomen and Pelvis
Sagittal T2 FS	26.0 cm	5 mm	1.5 mm	Pelvis
Axial T1 FS	25.0 cm	5 mm	1.0 mm	Pelvis
Axial T2 FS	25.0 cm	5 mm	1.0 mm	Pelvis
Check with Radiologist - may need orthogonal to vagina				
Coronal T2	24.0 cm	2.0 mm	0.5mm	Uterus/Cervix
Sagittal T2	24.0 cm	2.0 mm	0.5 mm	Uterus/Cervix
Axial T2	24.0 cm	2.0 mm	0.5 mm	To Uterus (uterus only)