

MRI Adnexal Mass 3T only

Application Examples: adnexal mass

Patient Prep: KY gel in vagina

FoV Coverage: Bifurcation to pubic symphysis

Sequence	FoV	Thickness	Gap	
Sagittal T2	24.0 cm	5 mm	1.5 mm	Hip to hip
Axial T1	26.0 cm	5 mm	1.5 mm	Entire Pelvis
Axial T2	26.0 cm	5.0 mm	1.5 mm	Entire pelvis
Axial In/Out phase	35.0 cm	4.8 mm	No Gap	
Axial Diffusion	28.0 cm	4.0 mm	No Gap	B0, 600, 800, 1000
Coronal T2 FS	24.0 cm	5.0 mm	1.5 mm	
Axial Vibe FS/Lava Flex	40.0 cm	4.4 mm	No Gap	
INJECT CONTRAST (1mmol/Kg)				
Axial Vibe/ Lava Flex	40.0 cm	4.4mm	No Gap	70 seconds after injection