

MRI Perianal Fistula 3T only

Application Examples: Perianal fistula

Sequence	FoV	Thickness	Gap	
Axial T1	34.0 cm	5.0 mm	1.0 mm	Aortic Bifurcation to pubic symphysis
Axial T2	34.0 cm	5.0 mm	1.0 mm	Aortic Bifurcation to pubic symphysis
Remainder of exam scan orthogonal to anal canal				
Sagittal T2	18.0 cm	3.0 mm	0.5 mm	Rectal area only
Coronal Stir/Coronal T2 FS	18.0 cm	3.0 mm	0.5 mm	Rectal area only
Axial Stir/Axial T2 FS	20.0 mm	4.0 mm	1.0 mm	Rectal area only
Axial Vibe FS/Lava Flex pre	40.0 mm	4.4 mm	No Gap	
Check images with Radiologist				
INJECT CONTRAST (1mmol/Kg) - Hand inject				
Axial Vibe FS/Lava Flex	40.0 cm	4.4 mm	No Gap	No delay
Coronal Vibe FS/Lava Flex	40.0 cm	4.4 mm	No Gap	