Gundersen Health System Ortho MR Protocol GE 1.5T 450W

Routine Elbow

Application Examples: Elbow pain, ulnar nerve entrapment, pitchers elbow

Sequence	FoV	Thickness	Gap
Axial PD	14.0	3mm	0.5mm
Axial T2 FS	14.0	3.0mm	0.5mm
Coronal T1	14.0	3.0mm	0.5mm
Coronal T2 FS	14.0	3.0mm	0.5mm
Sagittal T1	14.0	3.0mm	0.5mm
Sagittal T2 FS	14.0	3.0mm	0.5mm







Axial Coverage

Coronal Coverage

Sagittal Coverage